

## WHAT IS SILVRSTAT?

Silvrstat is a powerful combination of antibacterial wound care products and a healing skin cream. Silvrstat products are made using a patented technology referred to as Hydronanon Technology®. Hydronanon Technology is made up of metallic silver nanoparticles ( $\text{Ag}_0$ ) coated with tetrasilver tetroxide ( $\text{Ag}_4\text{O}_4$ ).



***Wound Care  
doesn't have to be  
complex, it can be  
SilvrSTAT simple.***

[info@ablmedical.com](mailto:info@ablmedical.com)

801-763-8000

[www.silvrstat.com](http://www.silvrstat.com)



## HYDRONANON TECHNOLOGY®

**Gentle 3-Step PTP Wound  
Care with SilvrSTAT  
Protocol**



PREPARE | TREAT | PROTECT

# INDICATIONS

## WOUND WASH

### For use of OTC Application for:

- NORMAL SKIN
- MINOR WOUNDS
- ABRADED SKIN
- IRRITATED AREAS

### For use under a licensed healthcare practitioner:

- STAGE I-IV PRESSURE ULCERS, STASIS ULCERS, DIABETIC ULCERS, FOOT ULCERS
- POST-SURGICAL WOUNDS
- 1ST AND 2ND DEGREE BURNS
- CUTS & ABRASIONS
- MINOR SKIN IRRITATIONS

## WOUND DRESSING GEL

### For use of OTC Application for:

- FIRST AND SECOND DEGREE BURNS
- MINOR CUTS
- SKIN IRRITATIONS
- LACERATIONS
- ABRASIONS

### For use under a licensed healthcare practitioner:

- WOUNDS SUCH AS STASIS ULCERS, PRESSURE ULCERS, DIABETIC ULCERS
- SKIN TEARS
- SURGICAL INCISION SITES
- DEVICE INSERTION SITE WOUNDS
- GRAFT SITES
- DONOR SITES


SILVRSTAT® ANTIBACTERIAL WOUND DRESSING GEL CONTAINS SILVER THAT IN LABORATORY TESTS HAS BEEN SHOWN TO INHIBIT THE GROWTH OF MICROORGANISMS SUCH AS: Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, antibiotic-resistant bacteria MRSA and VRE, as well as fungi such as Candida albicans.

# UNIQUE TECHNOLOGY


We have been serving doctors, physicians and patients for over ten years, with a mission of serving others with the power of our technology.

SilvrSTAT Hydronanon Technology®, was developed with the goal of providing a new generation of silver-based products that are safe and effective for use on most types of wounds.

Nano Silver Molecule



Ag<sub>4</sub>O<sub>4</sub> Molecular shell



• Oxygen  
• Silver (Ag)

# OUR PTP WOUND PROTOCOL

## A quick and gentle 3-step PTP process:

**Prepare-** Use the Nano Silver Wound Wash by spraying on wound.

**Treat-** Use the Wound Dressing Gel by rubbing a thin layer of gel directly over affected area.

**Protect-** Protect the area with the Skin Cream to make sure your hard work lasts!

